

OraQuick ADVANCE - Day of Test LOG

Test Site _____

Date _____

Site Type [] Clinic [] Outreach

Lab Technician _____

Test Kit LOT# _____

Test Kit EXP Date _____

Client Log

Client ID	Lab Temperature 15° - 37°C or 59° - 99°F	Start Time	Read Time	Was the test read 20-40 minutes after the start time? YES - test next CLT NO - retest this CLT	Result (✓ one)			Done by
					N	PP	I	
C Control HIV (-)				[] Yes [] No				
C Control HIV1 (+)				[] Yes [] No				
C Control HIV2 (+)				[] Yes [] No				
1				[] Yes [] No				
2				[] Yes [] No				
3				[] Yes [] No				
4				[] Yes [] No				
5				[] Yes [] No				
6				[] Yes [] No				
7				[] Yes [] No				
8				[] Yes [] No				
9				[] Yes [] No				
10				[] Yes [] No				

TOTAL - N

TOTAL - PP

TOTAL - I

RESULT = (N) Negative, (PP) Preliminary Positive, (I) Invalid

Notes _____

OraQuick ADVANCE