

OraQuick ADVANCE - Day of Test LOG

Test Site _____ Date _____
 Site Type []Clinic []Outreach Lab Technician _____
 Test Kit LOT# _____ Test Kit EXP Date _____

Client Log

Client ID	Lab Temperature 15° - 37°C or 59° - 99°F	Start Time	Read Time	Was the test read 20-40 minutes after the start time? YES - test next CLT NO - retest this CLT	Result (✓ one)			Done by
					N	PP	I	
C	Control HIV (-)			[]Yes []No				
C	Control HIV1 (+)			[]Yes []No				
C	Control HIV2 (+)			[]Yes []No				
1				[]Yes []No				
2				[]Yes []No				
3				[]Yes []No				
4				[]Yes []No				
5				[]Yes []No				
6				[]Yes []No				
7				[]Yes []No				
8				[]Yes []No				
9				[]Yes []No				
10				[]Yes []No				
11				[]Yes []No				
12				[]Yes []No				
13				[]Yes []No				
14				[]Yes []No				
15				[]Yes []No				
16				[]Yes []No				
17				[]Yes []No				
18				[]Yes []No				
19				[]Yes []No				
20				[]Yes []No				
21				[]Yes []No				
22				[]Yes []No				
23				[]Yes []No				
24				[]Yes []No				
25				[]Yes []No				
26				[]Yes []No				
27				[]Yes []No				
28				[]Yes []No				
29				[]Yes []No				
30				[]Yes []No				

TOTAL - N
 TOTAL - PP
 TOTAL - I

RESULT = (N) Negative, (PP) Preliminary Positive, (I) Invalid

Notes _____

O r a Q u i c k A D V A N C E